

Analysis and Feature Extraction from a PPG Signal

Purpose: The purpose of this project is to simulate, process, and analyze a synthetic PPG signal, which is commonly used in wearable devices to measure heart rate, respiratory rate, and heart rate variability (HRV) and other physiological parameters. The project demonstrates how to:

- Generate a synthetic PPG signal with noise.
- Clean and filter the signal to remove noise.
- Detect peaks and extract meaningful features from the signal.
- Perform time-domain and frequency-domain analysis
- Compute heart rate, respiratory rate, and heart rate variability (HRV).
- Explore signal processing techniques such as filtering, correlation, convolution, and frequency response analysis.

Worked Explanations:

1. Signal Generation and Preprocessing:

- A synthetic PPG signal is generated using a sinusoidal function with added noise to simulate real-world conditions
- The signal is visualized to understand its characteristics before and after adding noise.




Figure 01: Raw PPG Signal with noise and without noise

2. Noise Reduction and Filtering:

- A **bandpass filter** (0.5 – 3.0 Hz) is applied to remove unwanted noise and retain only the useful components of the PPG signal.
- A **low-pass filter** is also tested to observe its impact on the signal.




Figure 02: Noise reduced using bandpass filter




Figure 03: Noise reduced using bandpass filter

3. Signal Normalization:

- The PPG signal is normalized to a scale of 0 to 1 to enhance further analysis.




Figure 04: Normalized Signal

4. Peak Detection and Heart Rate Estimation:

- Peaks are detected using **NeuroKit2's** PPG find peaks function.
- The **valid peaks** are filtered based on amplitude thresholds.
- **Heart Rate (BPM)** is estimated using the detected peaks




Figure 05: Valid Peak Detect

5. Abnormality Detection:

- Inter-peak intervals (IBI) are analyzed.
- Abnormal peaks are identified based on deviations from the mean inter-beat interval.
- The percentage of abnormal peaks is calculated to assess irregularities.




Figure 06: Abnormal Peak

6. Feature Extraction:

- Statistical features of inter-beat intervals:
- **Mean, Standard Deviation (SDNN), Skewness, Kurtosis**
- **Abnormality percentage** in the detected peaks.
- **Frequency-domain analysis using FFT** (Fast Fourier Transform) to identify the dominant frequency in the signal.

Feature Extraction Results:

Mean Interval (s): **1.0025**

Standard Deviation of Intervals (s): **0.017139136501002624**

Skewness of Intervals: **0.3538001921438391**

Kurtosis of Intervals: **-1.4223630602082342**

Abnormality Percentage (%): **0.0**

Dominant Frequency (Hz): **1.0**

Heart Rate (BPM): **59.85037406483791**

Respiratory Rate (breaths per minute): **14.962593516209477)**

HRV Mean: **1.0025 seconds**

HRV SDNN: **0.017139136501002624 seconds**

Signal Energy: **123.67222088363499**

7. Signal Processing Techniques:

- **Correlation and convolution** are performed using a smoothing kernel to demonstrate signal processing techniques.
- **Cross-correlation** is computed to analyze the relationship between the original and delayed signals.
- **Auto-correlation** is computed to study the signal's self-similarity




Figure 07: Correlated Signal




Figure 08: Convolved Signal




Figure 09: Cross Correlation




Figure 10: Auto Cross Correlation

8. Filter Analysis:

- The impulse response, step response, phase response, and magnitude response of the bandpass filter are computed and visualized.




Figure 11: Impulse Response




Figure 12: Step Response




Figure 13: Phase Response




Figure 14: Magnitude Response

Outcomes and Applications:

- **Heart Rate (BPM) Estimation** – Helps in monitoring cardiovascular health.
- **Respiratory Rate Estimation** – Provides insights into breathing patterns.
- **Detection of Irregular Heartbeats** – Helps in identifying potential arrhythmias or heart abnormalities.
- **Feature Extraction for Machine Learning Models** – Can be used for automated health monitoring applications.
- **Signal Filtering Techniques Comparison** – Evaluates different noise reduction methods.
- **HRV Analysis** – Assesses autonomic nervous system activity and overall heart health.

Source Code:

```
import numpy as np
import neurokit2 as nk
import matplotlib.pyplot as plt
from scipy.signal import butter, filtfilt, correlate
from scipy.stats import skew, kurtosis, mode

# Generate the PPG Signal with Noise

# Generate a 10-second synthetic PPG signal with sinusoidal behavior and random noise.
fs = 100 # Sampling frequency (Hz)
t = np.linspace(0, 10, fs * 10) # Time vector for 10 seconds

ppg_signal = 1 + 0.5 * np.sin(2 * np.pi * 1 * t) + 0.1 * np.random.randn(len(t)) # PPG signal with noise

# Plot the raw PPG signal plt.figure(figsize=(10, 4))

plt.plot(t, ppg_signal, label="Raw PPG Signal")

plt.title("Raw PPG Signal with Noise")

plt.xlabel("Time (s)")

plt.ylabel("Amplitude")

plt.legend()

plt.show()
```

```

# Clean PPG signal (without noise)

ppg_signal_clean = 1 + 0.5 * np.sin(2 * np.pi * 1 * t) # PPG signal without noise

# Plot the clean PPG signal

plt.figure(figsize=(10, 4))

plt.plot(t, ppg_signal_clean, label="Raw PPG Signal", color='blue')

plt.title("Raw PPG Signal (Without Noise)")

plt.xlabel("Time (s)")

plt.ylabel("Amplitude")

plt.legend()

plt.grid()

plt.show()

from scipy.signal import butter, filtfilt

# Define bandpass filter

low_cutoff = 0.5 # Lower cutoff frequency (Hz)

high_cutoff = 3.0 # Upper cutoff frequency (Hz)

b, a = butter(4, [low_cutoff / (fs / 2), high_cutoff / (fs / 2)], btype='band')

# Apply the filter to the signal

filtered_signal = filtfilt(b, a, ppg_signal)

# Plot the filtered signal

plt.figure(figsize=(10, 4))

plt.subplot(121)

plt.plot(t, ppg_signal)

plt.title("Original signal ")

plt.xlabel("Time (s)")

plt.ylabel("Amplitude")

plt.subplot(122)

plt.plot(t, filtered_signal, color="orange")

```

```
plt.title("Noise-Reduced Using Bandpass Filter")  
plt.xlabel("Time (s)")  
plt.ylabel("Amplitude")  
plt.legend()  
plt.show()
```

```
# Define low-pass filter
```

```
cutoff = 3.0 # Cutoff frequency (Hz)
```

```
b, a = butter(4, cutoff / (fs / 2), btype='low')
```

```
# Apply the filter to the signal
```

```
low_passed_signal = filtfilt(b, a, ppg_signal)
```

```
# Plot the low-pass filtered signal
```

```
plt.figure(figsize=(10, 4))
```

```
plt.subplot(121)
```

```
plt.plot(t, ppg_signal)
```

```
plt.title("Orginal Signal ")
```

```
plt.xlabel("Time (s)")
```

```
plt.ylabel("Ampllidtue")
```

```
plt.subplot(122)
```

```
plt.plot(t, low_passed_signal, color="green")
```

```
plt.title("Noise-Reduced Using Low-pass Filter")
```

```
plt.xlabel("Time (s)")
```

```
plt.ylabel("Amplitude")
```

```
plt.legend()
```

```
plt.show()
```



```
# Normalize the PPG signal between 0 and 1.
```

```
normalized_signal = (ppg_signal - np.min(ppg_signal)) / (np.max(ppg_signal) -  
np.min(ppg_signal)) # Normalization
```

```
# Plot the normalized signal
```

```
plt.figure(figsize=(10, 4))
```

```
plt.plot(t, normalized_signal, label="Normalized Signal")
```

```
plt.title("Normalized Signal")
```

```
plt.xlabel("Time (s)")
```

```
plt.ylabel("Normalized Amplitude")
```

```
plt.legend()
```

```
plt.show()
```

```
# Detect peaks using nk.ppg_findpeaks
```

```
peaks = nk.ppg_findpeaks(filtered_signal, sampling_rate=fs)["PPG_Peaks"]
```

```
# Filter valid peaks based on amplitude thresholds
```

```
valid_peaks = peaks[(filtered_signal[peaks] > 0.2) & (filtered_signal[peaks] < 1.8)] # Valid peaks
```

```
# Plot the filtered signal with valid peaks
```

```
plt.figure(figsize=(10, 4))
```

```
plt.plot(t, filtered_signal, label="Filtered Signal", color='orange')
```

```
plt.scatter(t[valid_peaks], filtered_signal[valid_peaks], color='blue', label="Valid Peaks")
```

```
plt.title("Filtered PPG Signal with Valid Peaks")
```

```
plt.xlabel("Time (s)")
```

```
plt.ylabel("Amplitude")
```

```
plt.legend()
```

```
plt.show()
```

```
# Detect irregular peak intervals
```

```

inter_peak_intervals = np.diff(valid_peaks) / fs # Intervals in seconds

mean_ibi = np.mean(inter_peak_intervals)

std_ibi = np.std(inter_peak_intervals)

# Define thresholds for abnormal intervals (e.g., mean  $\pm$  2*std)

lower_threshold = mean_ibi - 2 * std_ibi

upper_threshold = mean_ibi + 2 * std_ibi

abnormal_intervals = (inter_peak_intervals < lower_threshold) | (inter_peak_intervals >
upper_threshold)

# Find indices of abnormal intervals

abnormal_peaks = valid_peaks[1:][abnormal_intervals]

# Plot detected abnormalities on the filtered signal

plt.figure(figsize=(10, 4))

plt.plot(t, filtered_signal, label="Filtered Signal", color='orange')

# plt.scatter(t[valid_peaks], filtered_signal[valid_peaks], color='green', label="Valid Peaks")

plt.scatter(t[abnormal_peaks], filtered_signal[abnormal_peaks], color='blue', label="Abnormal
Peaks")

plt.title("Filtered Signal with Abnormal Peaks")

plt.xlabel("Time (s)")

plt.ylabel("Amplitude")

plt.legend()

plt.show()

# Feature 1: Statistical features of inter-peak intervals

mean_interval = np.mean(inter_peak_intervals)

std_interval = np.std(inter_peak_intervals)

skewness_interval = skew(inter_peak_intervals)

```

```

kurtosis_interval = kurtosis(inter_peak_intervals)

# Feature 2: Abnormality percentage
abnormality_percentage = len(abnormal_peaks) / len(valid_peaks) * 100

# Feature 3: FFT-based frequency domain analysis
fft_signal = np.fft.fft(filtered_signal)
frequencies = np.fft.fftfreq(len(filtered_signal), 1 / fs)
dominant_frequency = frequencies[np.argmax(np.abs(fft_signal[:len(fft_signal) // 2]))]

# Print extracted features
print("Feature Extraction Results:")
print(f"Mean Interval (s): {mean_interval}")
print(f"Standard Deviation of Intervals (s): {std_interval}")
print(f"Skewness of Intervals: {skewness_interval}")
print(f"Kurtosis of Intervals: {kurtosis_interval}")
print(f"Abnormality Percentage (%): {abnormality_percentage}")
print(f"Dominant Frequency (Hz): {dominant_frequency}")

# Estimate heart rate and respiratory rate from valid peaks.
heart_rate = 60 / np.mean(np.diff(valid_peaks) / fs) # Calculate heart rate (BPM)
respiratory_rate = heart_rate / 4 # Approximate respiratory rate assuming 4:1 HR:RR ratio

# Print heart rate and respiratory rate
print(f"Heart Rate (BPM): {heart_rate}")
print(f"Respiratory Rate (breaths per minute): {respiratory_rate}")

# Perform correlation using the kernel on the PPG signal.
kernel = np.ones(5) / 5 # Smoothing kernel
correlated_signal = np.correlate(ppg_signal, kernel, mode='same') # Correlation

# Plot the correlated signal

```

```
plt.figure(figsize=(10, 4))  
plt.plot(t, correlated_signal, label="Correlated Signal")  
plt.title("Correlated Signal")  
plt.xlabel("Time (s)")  
plt.ylabel("Amplitude")  
plt.legend()  
plt.show()
```

```
# Perform convolution on the PPG signal using a smoothing kernel.  
convolved_signal = np.convolve(ppg_signal, kernel, mode='same') # Convolution  
# Plot the convolved signal  
plt.figure(figsize=(10, 4))  
plt.plot(t, ppg_signal, label="Original Signal")  
plt.plot(t, convolved_signal, label="Convolved Signal")  
plt.title("Convolved Signal")  
plt.xlabel("Time (s)")  
plt.ylabel("Amplitude")  
plt.legend()  
plt.show()  
plt.show()
```

```
# Compute the cross-correlation between the original and delayed signals.  
cross_corr = correlate(ppg_signal, delayed_signal, mode='full') # Cross-correlation  
lags = np.arange(-len(ppg_signal) + 1, len(ppg_signal)) # Lags for correlation  
# Plot the cross-correlation  
plt.figure(figsize=(10, 4))  
plt.plot(lags, cross_corr, label="Cross-correlation")  
plt.title("Cross-correlation")
```

```
plt.xlabel("Lag")
```

```
plt.ylabel("Correlation")
```

```
plt.legend()
```

```
plt.show()
```

```
# Auto-correlation
```

```
# Compute the auto-correlation of the PPG signal.
```

```
auto_corr = correlate(ppg_signal, ppg_signal, mode='full') # Auto-correlation
```

```
lags = np.arange(-len(ppg_signal) + 1, len(ppg_signal)) # Lags for correlation
```

```
# Plot the auto-correlation
```

```
plt.figure(figsize=(10, 4))
```

```
plt.plot(lags, auto_corr, label="Auto-correlation")
```

```
plt.title("Auto-correlation")
```

```
plt.xlabel("Lag")
```

```
plt.ylabel("Correlation")
```

```
plt.legend()
```

```
plt.show()
```

```
# Design a bandpass filter and compute its impulse response.
```

```
impulse = np.zeros_like(ppg_signal) # Create an impulse signal
```

```
impulse[len(impulse) // 2] = 1 # Set the center to 1
```

```
b, a = butter(4, [0.5 / (fs / 2), 3.0 / (fs / 2)], btype='band') # Bandpass filter design
```

```
impulse_response = filtfilt(b, a, impulse) # Compute impulse response
```

```
# Plot the impulse response
```

```
plt.figure(figsize=(10, 4))
```

```
plt.plot(impulse_response, label="Impulse Response")
```

```
plt.title("Impulse Response")
```

```
plt.xlabel("Samples")
```

```
plt.ylabel("Amplitude")
```

```
plt.legend()
```

```
plt.show()
```

```
# Step Response
```

```
# Compute the step response of the bandpass filter.
```

```
step = np.ones_like(ppg_signal) # Create a step signal
```

```
step_response = filtfilt(b, a, step) # Compute step response
```

```
# Plot the step response
```

```
plt.figure(figsize=(10, 4))
```

```
plt.plot(step_response, label="Step Response")
```

```
plt.title("Step Response")
```

```
plt.xlabel("Samples")
```

```
plt.ylabel("Amplitude")
```

```
plt.legend()
```

```
plt.show()
```

```
# Phase Response
```

```
# Compute and plot the phase response of the filter.
```

```
from scipy.signal import freqz
```

```
w, h = freqz(b, a, worN=8000) # Frequency response
```

```
plt.figure(figsize=(10, 4))
```

```
plt.plot(w / np.pi, np.angle(h), label="Phase Response")
```

```
plt.title("Phase Response")
```

```
plt.xlabel("Normalized Frequency ( $\pi$  rad/sample)")
```

```
plt.ylabel("Phase (radians)")
```

```
plt.legend()
```

```
plt.show()
```

```
# Magnitude Response
```

```
# Compute and plot the magnitude response of the filter.
```

```
plt.figure(figsize=(10, 4))
```

```
plt.plot(w / np.pi, np.abs(h), label="Magnitude Response")
```

```
plt.title("Magnitude Response")
```

```
plt.xlabel("Normalized Frequency ( $\pi$  rad/sample)")
```

```
plt.ylabel("Magnitude")
```

```
plt.legend()
```

```
plt.show()
```

```
rr_intervals = np.diff(valid_peaks) / fs # RR intervals in seconds
```

```
hrv_mean = np.mean(rr_intervals)
```

```
hrv_sdnn = np.std(rr_intervals) # Standard deviation of RR intervals
```

```
print(f"HRV Mean: {hrv_mean} seconds")
```

```
print(f"HRV SDNN: {hrv_sdnn} seconds")
```

```
signal_energy = np.sum(filtered_signal**2)
```

```
print(f"Signal Energy: {signal_energy}")
```